



URGENT MEDICAL DEVICE CORRECTION

GE Healthcare

Healthcare Systems
9900 Innovation Drive
Wauwatosa, WI 53226
USA

<Date of Letter Deployment>

GEHC Ref#32036

To: Manager Biomedical/Clinical Engineering
Nurse Manager, Labor & Delivery/ NICU
Manager, Respiratory Therapy/ Director of Risk Management

RE: Potential Safety Issue Involving GE Warmers with Nellcor SpO₂

GE Healthcare has recently become aware of a potential safety issue associated with the Nellcor SpO₂ alarm settings on the following GE Giraffe and Panda iRes Warmers:

- Giraffe Warmers integrated with Nellcor SpO₂
- Panda iRes Warmers integrated with Nellcor SpO₂
- Panda Freestanding Warmers integrated with Nellcor SpO₂
- Panda Wall-Mount Warmers integrated with Nellcor SpO₂
- Nellcor SpO₂ Upgrade kit

Please ensure that all potential users in your facility are made aware of this safety notification and the recommended actions.

Safety Issue

There may be a loss of internal communication that manages the user-adjusted Nellcor SpO₂ alarm settings of the affected Giraffe and Panda iRes Warmer systems with Nellcor SpO₂ listed above. If the issue exists, the alarms may not activate as expected, which can result in false positive or false negative saturation and pulse rate alarm notifications to the caregiver.

The patient's oxygen saturation and pulse rate values are accurate, as displayed; however, the displayed user-set alarm limits may be different than the actual alarm limits used for alarm activation. All other clinical functionality of the warmer is unaffected.

This alarm failure is caused by a software issue that does not detect a loss of internal communication that manages the user-adjusted Nellcor SpO₂ alarm settings.

Safety Instructions

You may continue to use your system provided you follow the GE recommended actions:

- 1) If in use with a patient, turn off the Nellcor SpO₂ function in the control panel and switch to an approved alternate form of SpO₂ measurement.
- 2) Prior to using the affected GE warmer with the Nellcor SpO₂ function, or resuming the use of this function with any patient, contact your hospital Bio-Med or qualified technician to check the Nellcor SpO₂ alarm functionality using either an adult Nellcor SpO₂ sensor or a Nellcor SpO₂ simulator following the applicable instructions in the Appendix (see pages 3-4).
- 3) If the test activates a 'Low Pulse Rate' alarm, the user-adjusted Nellcor SpO₂ alarm limits are functioning properly and the warmer can be placed back in clinical service.

PLEASE NOTE: Alarm functionality must be re-checked following every service event requiring access to the internal electronics and during annual preventative maintenance, until revised software is installed.

4) If the test does NOT activate a 'Low Pulse Rate' alarm, remove the warmer from clinical service and contact your local GE Healthcare Service Representative.

**Affected
Product
Details**

All Giraffe Warmers integrated with Nellcor SpO₂
Panda iRes Warmers integrated with Nellcor SpO₂
Panda Freestanding Warmers integrated with Nellcor SpO₂
Panda Wall-Mount Warmers integrated with Nellcor SpO₂
Nellcor SpO₂ Upgrade kit

**Product
Correction**

GE Healthcare will correct all affected systems with a software revision at no cost to you. A GE Healthcare service representative will contact you to arrange for this correction.

**Contact
Information**

If you have any questions or concerns regarding this notification, please call the following phone numbers:

Saudi Arabia Toll Free number: 8004292222

Saudi Arabia Service Center: SaudiArabiaServiceCenter@ge.com


You can also contact your local GE Healthcare Service Representative.

Please be assured that maintaining a high level of safety and quality is our highest priority. If you have any questions, please contact us immediately per the contact information above.

Sincerely,



James Dennison
Vice President QARA
GE Healthcare Systems



Douglas M. Hansell, M.D., MPH
Chief Medical Officer
GE Healthcare

Appendix

A. Testing Nellcor SpO₂ Alarm Functionality Using SpO₂ Simulator

To confirm that the warmer with Nellcor SpO₂ is activating the alarms based on the user-set limits, ensure patient is removed, and then generate a 'Low Pulse Rate' alarm by completing the following steps:

1. Connect SpO₂ simulator to the Nellcor probe jack shown below.



2. Turn On the SpO₂ feature by pressing the key to the left of the SpO₂ box.



3. Access the SpO₂ Alarm Limit settings by first pressing "Menu" key on the right of the screen and then the "SpO₂" key. Next, press the "Alarm Limits" key (shown below) to change the alarm limit settings for pulse rate or SpO₂.



4. Set the 'Low Pulse Rate' limit to 110.
5. Using the SpO₂ Simulator, simulate a pulse rate greater than 40 and lower than 110.
6. Confirm a "Low Pulse Rate" alarm activates.

Appendix

B. Testing Nellcor SpO₂ Alarm Functionality Using an Adult Nellcor SpO₂ Sensor

To confirm that the warmer with Nellcor SpO₂ is activating the alarms based on the user-set limits, ensure patient is removed, and then generate a 'Low Pulse Rate' alarm by completing the following steps:


1. Position the adult SpO₂ sensor on the adult finger.
2. Next, connect the sensor to the patient cable.
3. Plug the patient cable into the Nellcor probe jack shown below.



4. Turn On the SpO₂ feature by pressing the button to the left of the SpO₂ box.



5. Access the SpO₂ Alarm Limit settings by first pressing "Menu" key on the right of the screen and then the "SpO₂" key. Next, press the "Alarm Limits" key (shown below) to change the alarm limit settings for pulse rate or SpO₂.



6. Set the Low Pulse Rate limit to 110.
7. Measure the pulse rate. The measured pulse rate must be greater than 40 and less than 110, if it is not please use a SpO₂ simulator to test the system.
8. Confirm a **Low Pulse Rate** alarm activates.