

**To the ATTENTION of:**  
**Doctors using Synthes Mandible**  
**Distractor multiaxial**

01 October 2012

## URGENT: MEDICAL DEVICE PRODUCT REMOVAL

Part Description / Part Number:

Part Number	Part Description	Lot Number
487.931	Distractor Body, Titanium Alloy (TAV)	7895264
		7793390
		7583381
		7583380
		7528840
		3708643
		3683421
		3683403
487.966	Mandible Distractor, multiaxial, right, Titanium Alloy (TAV)	3560888
		7765627
		7562081
		7532790
		3800757
487.967	Mandible Distractor, multiaxial, left, Titanium Alloy (TAV)	3757207
		7834374
		7662788
		7562083
		3757213
		3739743

A



Dear Doctor

Synthes is initiating a Medical Device Product Removal related to the Synthes Mandible Distractor, multitaxial.

This action is being initiated following a detailed investigation in response to a reported complaint of missing outer sheer pin. The inner diameter (for fixing the sheer pin) of the worm gear component was out of specification.


For patients who may be undergoing treatment with this device, it is advised to contact the patient to check if the device they have is one of the effected lots. If the patient is in the consolidation phase of treatment, a consolidation rod should be used according to the Surgical Technique, page 13.

If the patient is in the distraction phase of treatment and has one of these device lots, please follow the steps on page 13 of the Surgical Technique Guide which describes the application of a consolidation rod and removal of the device. Once the central body of the device has been removed as described, a new device can be put in its place and secured. The distraction arms can be preloaded by normal activation. The consolidation rod can then be removed.

To eliminate any risks Synthes decided to remove all potential affected parts from the field. Synthes is requesting that you immediately cease using the product and please examine your inventory for the above effected lot numbers and remove them. When in doubt, please remove the product and provide it to your Synthes Sales Representative.

Thank you for your attention and cooperation.

Synthes GmbH

A handwritten signature in blue ink, appearing to be 'H. Gribi'.

Heinz Gribi  
Manager Complaint Handling Unit

A handwritten signature in blue ink, appearing to be 'M. Wien'.

Markus Wien  
Director Quality EMEA

Cc:

A handwritten mark in blue ink, possibly a stylized 'A' or 'F'.



**NOTICE: MEDICAL DEVICE RECALL****Mandible Distractor, multiaxial**  
**Verification Section**

Part Number	Part Description	Lot Number
487.931	Distractor Body, Titanium Alloy (TAV)	7895264
		7793390
		7583381
		7583380
		7528840
		3708643
		3683421
		3683403
487.966	Mandible Distractor, multiaxial, right, Titanium Alloy (TAV)	3560888
		7765627
		7562081
		7532790
		3800757
487.967	Mandible Distractor, multiaxial, left, Titanium Alloy (TAV)	3757207
		7834374
		7662788
		7562083
		3757213
		3739743

☐ We have located the identified product in stock; returned quantity is documented below, and have retained a copy of this letter for our records.

☐ We do not have any identified product in stock; returned quantity is zero, and have retained a copy of this letter for our records.

RETURNED DEVICES (including quantity):

---

---

---

---

Name/Title (please print) \_\_\_\_\_

Phone Number: \_\_\_\_\_

Signature and Date: \_\_\_\_\_

