



Mortara Instrument, Inc.

7865 North 86th Street
Milwaukee, WI 53224

Tel: 414.354.1600

Fax: 414.354.4760

www.mortara.com

URGENT MEDICAL DEVICE CORRECTION

October 3, 2014

To: Health Care Administrator/ Risk Manager
Director of Biomedical/ Clinical Engineering

Re: Improper mounting of Mortara Surveyor S12 and Surveyor S19 Patient Monitors may result in a safety hazard.

Dear Customer,

Mortara Instrument, Inc. has recently become aware of a potential safety hazard involving our Surveyor S12 and Surveyor S19 Patient Monitors. A hazardous situation may be created when mounting these monitors using the 75mm VESA mounting hole pattern provided in the battery cover of these monitors. If the mounting screws provided by Mortara are not used and if the mounting screws used are longer than the provided screws, it is then possible for the longer screws to be driven completely through the battery cover and then possibly penetrate the protective casing of the lithium-ion battery which may result in a fire.

There have been no reported incidents of fire, injury or death due to this potential safety hazard.

A correction has been developed to eliminate the potential risk. The correction replaces the original battery cover with a new design incorporating a barrier that prevents mounting screws from penetrating into the battery compartment.

Please ensure that all potential users in your facility are made aware of this safety notification and that the required actions listed below are taken immediately.

1. Locate each Surveyor S12 and Surveyor S19 unit. The list of units delivered by Mortara to your organization is included.
2. For each unit, follow the included instructions (document number MIS-11-183-01) detailing the correction steps.
3. Complete and return the attached Medical Device Correction Return Response Acknowledgement and Receipt Form. A response is required for each device listed on the form.



Please be assured that maintaining a high level of safety and quality is our highest priority. If you have any questions regarding this notice, please contact Technical Support at +1.414.354.1600 or TechSupport@mortara.com , or contact your local Service Representative.

Sincerely,

Timothy E. Field
Senior Vice President Quality



MEDICAL DEVICE CORRECTION RETURN RESPONSE
Acknowledgement and Receipt Form
 Response is Required

Customer Information:

**SAUDI HEALTH SERVICE CO,LTD
 TASLEEH STREET, NEAR GIRLS SCHOOL
 GUAIZA DISTRICT
 JEDDAH, WR .. 21443
 SAUDI ARABIA**

Surveyor S12 and Surveyor S19 Patient Monitor Medical Device Correction

I have read and understand the recall instructions provided with the October 3, 2014 letter. Yes No

Were any adverse events associated with the recalled product? Yes No

If yes, please explain:

Affected Product Information by Serial number:

Product Name	Mortara Part Number	Device Serial Number	Correction Completed
SURVEYOR	SUR12-BCF-AXXAD	114170207572	YES NO
SURVEYOR	SUR19-BCF-AXXAC	114170207577	YES NO

Signature of Receipt and Completion _____

Name/Title	
Telephone	
Email address	

FAX completed form to: +1.414.354.4760, ATTN: S12/S19 CORRECTION

OR e-mail completed form to: CORRECTIONS@MORTARA.COM

OR mail completed form to:

MORTARA INSTRUMENT, INC.
 ATTN: S12/S19 CORRECTION
 7865 NORTH 86TH STREET
 MILWAUKEE, WI 53224