


Field Safety Notice

Urgent Medical Device Correction - 2955842-02/01/16-003-C

da Vinci® Si™ and Xi™ Touchpad Compact Flash Card

<p>Introduction and Reason for Field Action</p>	<p>Dear <i>da Vinci</i> Customer,</p> <p>This Field Safety Notice is to advise you that Intuitive Surgical is initiating a voluntary correction to the <i>da Vinci</i> Si and Xi Surgical System, relating to a component used in the Surgeon Console and Patient Side Carts.</p> <p>The surgeon side console (SSC) of the Si and Xi systems and the patient side cart (PSC) of the Xi systems have touchpads that contain compact flash cards that are used to save user preferences and to assist in system setup (Figure 1 below).</p> <p>Intuitive Surgical has found that compact flash cards used in touchpads on specific <i>da Vinci</i> Si and Xi Surgical Systems have the potential to become corrupted when the touchpads experience an abrupt power loss and/or are powered on and off within a short period of time.</p> <p>The corruption of the compact flash occurs when power is cycled. If the compact flash has been corrupted, the <i>da Vinci</i> system will display a fault during the startup of the system. In some cases the faults may be recoverable, which allow continued use of the system with the touchpad functionality disabled. In other cases the faults may be non-recoverable and will prevent use of the system.</p> <p>The root cause for this issue was a change made by the manufacturer of the compact flash card. After discovery of the issue, a new, more robust design was implemented in November 2015. All Surgeon Consoles and Patient Side carts manufactured after November are not affected by this issue.</p>  <p>Figure 1. Touchpad locations on the <i>da Vinci</i> Xi Patient Side Cart (Left) and <i>da Vinci</i> Si and Xi Surgeon Console (Right).</p>
<p>Risk to Health</p>	<p>Through December 2015, 8.1% of affected touchpads have failed.</p> <p>There have been no reported injuries or deaths as a result of this issue. In the worst case scenario, the <i>da Vinci</i> surgical system will encounter a non-recoverable fault, preventing use of the system. The error typically presents itself during startup, so will</p>

	not likely occur during a procedure (unless a mid-procedure power restart is required). If it is required to convert a procedure to an alternate surgical method due to the fault, the patient may be at higher risk of surgical complications because of the historically higher rate of complications in open surgery.									
<p>Affected Countries and Products</p>	<p>Affected Countries: Australia, Austria, Belgium, Brazil, Canada, Chile, China, Czech Republic, Denmark, Finland, France, Germany, Great Britain, Hong Kong, India, Israel, Italy, Japan, Malaysia, Monaco, Netherlands, Philippines, Portugal, Qatar, Romania, Russia, Saudi Arabia, South Korea, Sweden, Switzerland, Taiwan, Turkey, United States, and Venezuela.</p> <table border="1" data-bbox="461 620 1131 721"> <thead> <tr> <th data-bbox="461 620 1131 654">Affected Product</th> </tr> </thead> <tbody> <tr> <td data-bbox="461 654 1131 687"><i>da Vinci Xi Surgical System</i></td> </tr> <tr> <td data-bbox="461 687 1131 721"><i>da Vinci Si Surgical System</i></td> </tr> </tbody> </table> <table border="1" data-bbox="461 763 1131 898"> <thead> <tr> <th data-bbox="461 763 746 831">Affected Component Part Number(s)</th> <th data-bbox="754 763 1131 831">Part Name(s)</th> </tr> </thead> <tbody> <tr> <td data-bbox="461 831 746 864">187505-04</td> <td data-bbox="754 831 1131 864"><i>da Vinci Si/Xi SSC Touchpad</i></td> </tr> <tr> <td data-bbox="461 864 746 898">187716-03</td> <td data-bbox="754 864 1131 898"><i>da Vinci Xi PSC Touchpad</i></td> </tr> </tbody> </table>	Affected Product	<i>da Vinci Xi Surgical System</i>	<i>da Vinci Si Surgical System</i>	Affected Component Part Number(s)	Part Name(s)	187505-04	<i>da Vinci Si/Xi SSC Touchpad</i>	187716-03	<i>da Vinci Xi PSC Touchpad</i>
Affected Product										
<i>da Vinci Xi Surgical System</i>										
<i>da Vinci Si Surgical System</i>										
Affected Component Part Number(s)	Part Name(s)									
187505-04	<i>da Vinci Si/Xi SSC Touchpad</i>									
187716-03	<i>da Vinci Xi PSC Touchpad</i>									
<p>Actions to be taken by the Customer/ User</p>	<ul style="list-style-type: none"> • If your systems display any faults, follow the standard troubleshooting procedure. • Distribute a copy of this letter to all <i>da Vinci</i> users at your facility. 									
<p>Actions to be taken by Intuitive Surgical</p>	<p>An Intuitive Surgical representative will schedule a site visit to replace the affected compact flash card at your site.</p> <p>Note that replacing the compact flash card may require users to re-enter their user preferences.</p>									
<p>Further Information & Support</p>	<p>If you need further information or support concerning this Medical Device Notification, please contact your Clinical Sales Representative or contact Intuitive Surgical Customer Service at the numbers listed below:</p> <ul style="list-style-type: none"> • North and South America: (800) 876-1310, Option 3 (6 AM to 5 PM PST) or mail: customersupport-servicesupport@intusurg.com • Europe, Middle East, Asia and Africa: +800 0821 2020 or +41 21 821 2020 (8 AM to 6 PM CET) or ics@intusurg.com • South Korea: 02-3271-3200 (9 AM to 6 PM KSTJ) • Japan: 0120-56-5635 or 03-5575-1362 (9 AM to 6 PM JST) 									

Please be informed that the appropriate Regulatory Authority for your region has been notified of this Medical Device Correction.

Sincerely,

Intuitive Surgical Sàrl

Chemin des Mûriers 1

CH-1170 Aubonne, Switzerland

+41 21 821 2020