

URGENT NOTICE: MEDICAL DEVICE RECALL
Synthes Trochanteric Fixation Nailing (TFN) System

Dear Sir/ Madam:

Synthes is initiating a voluntary recall of all lots of the Trochanteric Fixation Nailing (TFN) System due to the possibility that the outer pouch may be compromised. This recall is not being undertaken on the basis of adverse events.

The issue was identified during a validation test for a change to our sterile packaging. The test revealed there is the potential for damage to the external pouch. The product is packaged in two pouches; a sterile, internal pouch within an external pouch. While the risk of harm to the patient is remote, we are recalling the products as a measure of caution. An improvement to the packaging line has been implemented, and is now in place to help prevent the issue from arising again.

Synthes is requesting that you please examine your inventory for all products with the Part numbers listed below. Synthes strongly urges you to return the identified devices and replacement product will be provided once it becomes available.

If you have any questions or concerns, please call your Synthes Sales organization.

November 1st, 2012

Return Receipt Requested

SYNTHES GmbH
Luzernstrasse 21
CH-4528 Zuchwil
Tel +41 61 965 11 11

FULL RECALLED PRODUCT LIST:

Part Description	Part Number	Lot Number
Trochanteric Fixation Nails	See Attached	See Attached-Please note, All Lot numbers up to and including Lot # 7060897 are within the scope of this recall.

If you **DO have** any of the identified devices and intend to return them for a replacement, please take the following steps:

- Complete the Verification Section at the end of this letter by checking the appropriate box indicating affected product has been located. Also, please indicate the number of devices found. Please include your name, title, telephone number and signature in the spaces provided.
- Return the Verification Form (page 3 of this letter) with the product to your local Synthes Sales organization.

If you **DO NOT** have the identified product, please take the following steps:

- Complete the attached Verification Section at the end of this letter by checking the appropriate box indicating that no affected product has been located. Please include your name, title, telephone number and signature in the spaces provided. This return documentation acknowledges your receipt of medical device removal information.
- Return the document to your Synthes Sales organization.

If you have any questions, please contact your Synthes Sales Representative.

We apologize for any inconvenience this may cause and thank you in advance for your attention and cooperation. Please contact your Synthes Sales Consultant for additional information as required.

Sincerely,



Heinz Gribi
Manager, Complaint Handling Unit
Synthes GmbH

URGENT NOTICE: MEDICAL DEVICE RECALL
Synthes Trochanteric Fixation Nailing (TFN) System

Verification Section

Part Description

Part Description	Part Number	Lot Number
Trochanteric Fixation Nails	See Attached	See Attached-Please note, All Lot numbers up to and including Lot # 7060897 are within the scope of this recall

- I have located the identified product in stock; returned quantity is documented below for replacement.
- I do not have any identified product in stock; returned quantity is zero.

RETURNED DEVICES (including quantity):

Name/Title: (please print)

Phone Number:

Signature and Date:

Product Number	Product Description	Lot Number
456.314S	10MM/125 DEG TI CANN TROCH FIXATION NAIL 170MM-STERILE	All
456.315S	10MM/130 DEG TI CANN TROCH FIXATION NAIL 170MM-STERILE	All
456.316S	10MM/135 DEG TI CANN TROCH FIXATION NAIL 170MM-STERILE	All
456.317S	11MM/125 DEG TI CANN TROCH FIXATION NAIL 170MM-STERILE	All
456.318S	11MM/130 DEG TI CANN TROCH FIXATION NAIL 170MM-STERILE	All
456.319S	11MM/135 DEG TI CANN TROCH FIXATION NAIL 170MM-STERILE	All
456.321S	12MM/125 DEG TI CANN TROCH FIXATION NAIL 170MM-STERILE	All
456.322S	12MM/130 DEG TI CANN TROCH FIXATION NAIL 170MM-STERILE	All
456.323S	12MM/135 DEG TI CANN TROCH FIXATION NAIL 170MM-STERILE	All
456.324S	10MM/125 DEG TI CANN TROCH FIXATION NAIL 235MM-STERILE	All
456.325S	10MM/130 DEG TI CANN TROCH FIXATION NAIL 235MM-STERILE	All
456.326S	10MM/135 DEG TI CANN TROCH FIXATION NAIL 235MM-STERILE	All
456.327S	11MM/125 DEG TI CANN TROCH FIXATION NAIL 235MM-STERILE	All
456.328S	11MM/130 DEG TI CANN TROCH FIXATION NAIL 235MM-STERILE	All
456.329S	11MM/135 DEG TI CANN TROCH FIXATION NAIL 235MM-STERILE	All
456.330S	10MM/125 DEG TI CANN TROCH FIXATION NAIL 300MM/RIGHT-STER	All
456.331S	10MM/125 DEG TI CANN TROCH FIXATION NAIL 300MM/LEFT-STER	All
456.332S	10MM/125 DEG TI CANN TROCH FIXATION NAIL 320MM/RIGHT-STER	All
456.333S	10MM/125 DEG TI CANN TROCH FIXATION NAIL 320MM/LEFT-STER	All
456.334S	10MM/125 DEG TI CANN TROCH FIXATION NAIL 340MM/RIGHT-STER	All
456.335S	10MM/125 DEG TI CANN TROCH FIXATION NAIL 340MM/LEFT-STER	All
456.336S	10MM/125 DEG TI CANN TROCH FIXATION NAIL 360MM/RIGHT-STER	All
456.337S	10MM/125 DEG TI CANN TROCH FIXATION NAIL 360MM/LEFT-STER	All
456.338S	10MM/125 DEG TI CANN TROCH FIXATION NAIL 380MM/RIGHT-STER	All
456.339S	10MM/125 DEG TI CANN TROCH FIXATION NAIL 380MM/LEFT-STER	All
456.340S	10MM/125 DEG TI CANN TROCH FIXATION NAIL 400MM/RIGHT-STER	All
456.341S	10MM/125 DEG TI CANN TROCH FIXATION NAIL 400MM/LEFT-STER	All
456.342S	10MM/125 DEG TI CANN TROCH FIXATION NAIL 420MM/RIGHT-STER	All
456.343S	10MM/125 DEG TI CANN TROCH FIXATION NAIL 420MM/LEFT-STER	All
456.345S	10MM/125 DEG TI CANN TROCH FIXATION NAIL 440MM/LEFT-STER	All
456.346S	10MM/125 DEG TI CANN TROCH FIXATION NAIL 460MM/RIGHT-STER	All
456.347S	10MM/125 DEG TI CANN TROCH FIXATION NAIL 460MM/LEFT-STER	All
456.350S	10MM/130 DEG TI CANN TROCH FIXATION NAIL 300MM/RIGHT-STER	All
456.351S	10MM/130 DEG TI CANN TROCH FIXATION NAIL 300MM/LEFT-STER	All
456.352S	10MM/130 DEG TI CANN TROCH FIXATION NAIL 320MM/RIGHT-STER	All
456.353S	10MM/130 DEG TI CANN TROCH FIXATION NAIL 320MM/LEFT-STER	All
456.354S	10MM/130 DEG TI CANN TROCH FIXATION NAIL 340MM/RIGHT-STER	All
456.355S	10MM/130 DEG TI CANN TROCH FIXATION NAIL 340MM/LEFT-STER	All

456.356S	10MM/130 DEG TI CANN TROCH FIXATION NAIL 360MM/RIGHT-STER	All
456.357S	10MM/130 DEG TI CANN TROCH FIXATION NAIL 360MM/LEFT-STER	All
456.358S	10MM/130 DEG TI CANN TROCH FIXATION NAIL 380MM/RIGHT-STER	All
456.359S	10MM/130 DEG TI CANN TROCH FIXATION NAIL 380MM/LEFT-STER	All
456.360S	10MM/130 DEG TI CANN TROCH FIXATION NAIL 400MM/RIGHT-STER	All
456.361S	10MM/130 DEG TI CANN TROCH FIXATION NAIL 400MM/LEFT-STER	All
456.362S	10MM/130 DEG TI CANN TROCH FIXATION NAIL 420MM/RIGHT-STER	All
456.363S	10MM/130 DEG TI CANN TROCH FIXATION NAIL 420MM/LEFT-STER	All
456.364S	10MM/130 DEG TI CANN TROCH FIXATION NAIL 440MM/RIGHT-STER	All
456.365S	10MM/130 DEG TI CANN TROCH FIXATION NAIL 440MM/LEFT-STER	All
456.366S	10MM/130 DEG TI CANN TROCH FIXATION NAIL 460MM/RIGHT-STER	All
456.367S	10MM/130 DEG TI CANN TROCH FIXATION NAIL 460MM/LEFT-STER	All
456.370S	10MM/135 DEG TI CANN TROCH FIXATION NAIL 300MM/RIGHT-STER	All
456.371S	10MM/135 DEG TI CANN TROCH FIXATION NAIL 300MM/LEFT-STER	All
456.372S	10MM/135 DEG TI CANN TROCH FIXATION NAIL 320MM/RIGHT-STER	All
456.373S	10MM/135 DEG TI CANN TROCH FIXATION NAIL 320MM/LEFT-STER	All
456.374S	10MM/135 DEG TI CANN TROCH FIXATION NAIL 340MM/RIGHT-STER	All
456.375S	10MM/135 DEG TI CANN TROCH FIXATION NAIL 340MM/LEFT-STER	All
456.376S	10MM/135 DEG TI CANN TROCH FIXATION NAIL 360MM/RIGHT-STER	All
456.377S	10MM/135 DEG TI CANN TROCH FIXATION NAIL 360MM/LEFT-STER	All
456.378S	10MM/135 DEG TI CANN TROCH FIXATION NAIL 380MM/RIGHT-STER	All
456.379S	10MM/135 DEG TI CANN TROCH FIXATION NAIL 380MM/LEFT-STER	All
456.380S	10MM/135 DEG TI CANN TROCH FIXATION NAIL 400MM/RIGHT-STER	All
456.381S	10MM/135 DEG TI CANN TROCH FIXATION NAIL 400MM/LEFT-STER	All
456.382S	10MM/135 DEG TI CANN TROCH FIXATION NAIL 420MM/RIGHT-STER	All
456.383S	10MM/135 DEG TI CANN TROCH FIXATION NAIL 420MM/LEFT-STER	All
456.384S	10MM/135 DEG TI CANN TROCH FIXATION NAIL 440MM/RIGHT-STER	All
456.385S	10MM/135 DEG TI CANN TROCH FIXATION NAIL 440MM/LEFT-STER	All
456.386S	10MM/135 DEG TI CANN TROCH FIXATION NAIL 460MM/RIGHT-STER	All
456.387S	10MM/135 DEG TI CANN TROCH FIXATION NAIL 460MM/LEFT-STER	All
456.390S	11MM/125 DEG TI CANN TROCH FIXATION NAIL 300MM/RIGHT-STER	All
456.391S	11MM/125 DEG TI CANN TROCH FIXATION NAIL 300MM/LEFT-STER	All
456.392S	11MM/125 DEG TI CANN TROCH FIXATION NAIL 320MM/RIGHT-STER	All
456.393S	11MM/125 DEG TI CANN TROCH FIXATION NAIL 320MM/LEFT-STER	All
456.394S	11MM/125 DEG TI CANN TROCH FIXATION NAIL 340MM/RIGHT-STER	All
456.395S	11MM/125 DEG TI CANN TROCH FIXATION NAIL 340MM/LEFT-STER	All
456.396S	11MM/125 DEG TI CANN TROCH FIXATION NAIL 360MM/RIGHT-STER	All
456.397S	11MM/125 DEG TI CANN TROCH FIXATION NAIL 360MM/LEFT-STER	All

456.398S	11MM/125 DEG TI CANN TROCH FIXATION NAIL 380MM/RIGHT-STER	All
456.399S	11MM/125 DEG TI CANN TROCH FIXATION NAIL 380MM/LEFT-STER	All
456.400S	11MM/125 DEG TI CANN TROCH FIXATION NAIL 400MM/RIGHT-STER	All
456.401S	11MM/125 DEG TI CANN TROCH FIXATION NAIL 400MM/LEFT-STER	All
456.402S	11MM/125 DEG TI CANN TROCH FIXATION NAIL 420MM/RIGHT-STER	All
456.403S	11MM/125 DEG TI CANN TROCH FIXATION NAIL 420MM/LEFT-STER	All
456.404S	11MM/125 DEG TI CANN TROCH FIXATION NAIL 440MM/RIGHT-STER	All
456.405S	11MM/125 DEG TI CANN TROCH FIXATION NAIL 440MM/LEFT-STER	All
456.406S	11MM/125 DEG TI CANN TROCH FIXATION NAIL 460MM/RIGHT-STER	All
456.407S	11MM/125 DEG TI CANN TROCH FIXATION NAIL 460MM/LEFT-STER	All
456.410S	11MM/130 DEG TI CANN TROCH FIXATION NAIL 300MM/RIGHT-STER	All
456.411S	11MM/130 DEG TI CANN TROCH FIXATION NAIL 300MM/LEFT-STER	All
456.412S	11MM/130 DEG TI CANN TROCH FIXATION NAIL 320MM/RIGHT-STER	All
456.413S	11MM/130 DEG TI CANN TROCH FIXATION NAIL 320MM/LEFT-STER	All
456.414S	11MM/130 DEG TI CANN TROCH FIXATION NAIL 340MM/RIGHT-STER	All
456.415S	11MM/130 DEG TI CANN TROCH FIXATION NAIL 340MM/LEFT-STER	All
456.416S	11MM/130 DEG TI CANN TROCH FIXATION NAIL 360MM/RIGHT-STER	All
456.417S	11MM/130 DEG TI CANN TROCH FIXATION NAIL 360MM/LEFT-STER	All
456.418S	11MM/130 DEG TI CANN TROCH FIXATION NAIL 380MM/RIGHT-STER	All
456.419S	11MM/130 DEG TI CANN TROCH FIXATION NAIL 380MM/LEFT-STER	All
456.420S	11MM/130 DEG TI CANN TROCH FIXATION NAIL 400MM/RIGHT-STER	All
456.421S	11MM/130 DEG TI CANN TROCH FIXATION NAIL 400MM/LEFT-STER	All
456.422S	11MM/130 DEG TI CANN TROCH FIXATION NAIL 420MM/RIGHT-STER	All
456.423S	11MM/130 DEG TI CANN TROCH FIXATION NAIL 420MM/LEFT-STER	All
456.424S	11MM/130 DEG TI CANN TROCH FIXATION NAIL 440MM/RIGHT-STER	All
456.425S	11MM/130 DEG TI CANN TROCH FIXATION NAIL 440MM/LEFT-STER	All
456.426S	11MM/130 DEG TI CANN TROCH FIXATION NAIL 460MM/RIGHT-STER	All
456.427S	11MM/130 DEG TI CANN TROCH FIXATION NAIL 460MM/LEFT-STER	All
456.430S	11MM/135 DEG TI CANN TROCH FIXATION NAIL 300MM/RIGHT-STER	All
456.431S	11MM/135 DEG TI CANN TROCH FIXATION NAIL 300MM/LEFT-STER	All
456.432S	11MM/135 DEG TI CANN TROCH FIXATION NAIL 320MM/RIGHT-STER	All
456.433S	11MM/135 DEG TI CANN TROCH FIXATION NAIL 320MM/LEFT-STER	All
456.434S	11MM/135 DEG TI CANN TROCH FIXATION NAIL 340MM/RIGHT-STER	All
456.435S	11MM/135 DEG TI CANN TROCH FIXATION NAIL 340MM/LEFT-STER	All
456.436S	11MM/135 DEG TI CANN TROCH FIXATION NAIL 360MM/RIGHT-STER	All
456.437S	11MM/135 DEG TI CANN TROCH FIXATION NAIL 360MM/LEFT-STER	All
456.438S	11MM/135 DEG TI CANN TROCH FIXATION NAIL 380MM/RIGHT-STER	All
456.439S	11MM/135 DEG TI CANN TROCH FIXATION NAIL 380MM/LEFT-STER	All
456.440S	11MM/135 DEG TI CANN TROCH FIXATION NAIL 400MM/RIGHT-STER	All

456.441S	11MM/135 DEG TI CANN TROCH FIXATION NAIL 400MM/LEFT-STER	All
456.442S	11MM/135 DEG TI CANN TROCH FIXATION NAIL 420MM/RIGHT-STER	All
456.443S	11MM/135 DEG TI CANN TROCH FIXATION NAIL 420MM/LEFT-STER	All
456.444S	11MM/135 DEG TI CANN TROCH FIXATION NAIL 440MM/RIGHT-STER	All
456.445S	11MM/135 DEG TI CANN TROCH FIXATION NAIL 440MM/LEFT-STER	All
456.446S	11MM/135 DEG TI CANN TROCH FIXATION NAIL 460MM/RIGHT-STER	All
456.447S	11MM/135 DEG TI CANN TROCH FIXATION NAIL 460MM/LEFT-STER	All
456.450S	12MM/125 DEG TI CANN TROCH FIXATION NAIL 300MM/RIGHT-STER	All
456.451S	12MM/125 DEG TI CANN TROCH FIXATION NAIL 300MM/LEFT-STER	All
456.452S	12MM/125 DEG TI CANN TROCH FIXATION NAIL 320MM/RIGHT-STER	All
456.453S	12MM/125 DEG TI CANN TROCH FIXATION NAIL 320MM/LEFT-STER	All
456.454S	12MM/125 DEG TI CANN TROCH FIXATION NAIL 340MM/RIGHT-STER	All
456.455S	12MM/125 DEG TI CANN TROCH FIXATION NAIL 340MM/LEFT-STER	All
456.456S	12MM/125 DEG TI CANN TROCH FIXATION NAIL 360MM/RIGHT-STER	All
456.457S	12MM/125 DEG TI CANN TROCH FIXATION NAIL 360MM/LEFT-STER	All
456.458S	12MM/125 DEG TI CANN TROCH FIXATION NAIL 380MM/RIGHT-STER	All
456.459S	12MM/125 DEG TI CANN TROCH FIXATION NAIL 380MM/LEFT-STER	All
456.460S	12MM/125 DEG TI CANN TROCH FIXATION NAIL 400MM/RIGHT-STER	All
456.461S	12MM/125 DEG TI CANN TROCH FIXATION NAIL 400MM/LEFT-STER	All
456.462S	12MM/125 DEG TI CANN TROCH FIXATION NAIL 420MM/RIGHT-STER	All
456.463S	12MM/125 DEG TI CANN TROCH FIXATION NAIL 420MM/LEFT-STER	All
456.464S	12MM/125 DEG TI CANN TROCH FIXATION NAIL 440MM/RIGHT-STER	All
456.465S	12MM/125 DEG TI CANN TROCH FIXATION NAIL 440MM/LEFT-STER	All
456.466S	12MM/125 DEG TI CANN TROCH FIXATION NAIL 460MM/RIGHT-STER	All
456.467S	12MM/125 DEG TI CANN TROCH FIXATION NAIL 460MM/LEFT-STER	All
456.470S	12MM/130 DEG TI CANN TROCH FIXATION NAIL 300MM/RIGHT-STER	All
456.471S	12MM/130 DEG TI CANN TROCH FIXATION NAIL 300MM/LEFT-STER	All
456.472S	12MM/130 DEG TI CANN TROCH FIXATION NAIL 320MM/RIGHT-STER	All
456.473S	12MM/130 DEG TI CANN TROCH FIXATION NAIL 320MM/LEFT-STER	All
456.474S	12MM/130 DEG TI CANN TROCH FIXATION NAIL 340MM/RIGHT-STER	All
456.475S	12MM/130 DEG TI CANN TROCH FIXATION NAIL 340MM/LEFT-STER	All
456.476S	12MM/130 DEG TI CANN TROCH FIXATION NAIL 360MM/RIGHT-STER	All
456.477S	12MM/130 DEG TI CANN TROCH FIXATION NAIL 360MM/LEFT-STER	All
456.479S	12MM/130 DEG TI CANN TROCH FIXATION NAIL 380MM/LEFT-STER	All
456.480S	12MM/130 DEG TI CANN TROCH FIXATION NAIL 400MM/RIGHT-STER	All
456.481S	12MM/130 DEG TI CANN TROCH FIXATION NAIL 400MM/LEFT-STER	All
456.482S	12MM/130 DEG TI CANN TROCH FIXATION NAIL 420MM/RIGHT-STER	All
456.483S	12MM/130 DEG TI CANN TROCH FIXATION NAIL 420MM/LEFT-STER	All
456.484S	12MM/130 DEG TI CANN TROCH FIXATION NAIL 440MM/RIGHT-STER	All

November 1st, 2012

Return Receipt Requested

SYNTHES GmbH
Luzernstrasse 21
CH-4528 Zuchwil
Tel +41 61 965 11 11

456.485S	12MM/130 DEG TI CANN TROCH FIXATION NAIL 440MM/LEFT-STER	All
456.486S	12MM/130 DEG TI CANN TROCH FIXATION NAIL 460MM/RIGHT-STER	All
456.487S	12MM/130 DEG TI CANN TROCH FIXATION NAIL 460MM/LEFT-STER	All
456.490S	12MM/135 DEG TI CANN TROCH FIXATION NAIL 300MM/RIGHT-STER	All
456.491S	12MM/135 DEG TI CANN TROCH FIXATION NAIL 300MM/LEFT-STER	All
456.492S	12MM/135 DEG TI CANN TROCH FIXATION NAIL 320MM/RIGHT-STER	All
456.493S	12MM/135 DEG TI CANN TROCH FIXATION NAIL 320MM/LEFT-STER	All
456.494S	12MM/135 DEG TI CANN TROCH FIXATION NAIL 340MM/RIGHT-STER	All
456.495S	12MM/135 DEG TI CANN TROCH FIXATION NAIL 340MM/LEFT-STER	All
456.496S	12MM/135 DEG TI CANN TROCH FIXATION NAIL 360MM/RIGHT-STER	All
456.497S	12MM/135 DEG TI CANN TROCH FIXATION NAIL 360MM/LEFT-STER	All
456.498S	12MM/135 DEG TI CANN TROCH FIXATION NAIL 380MM/RIGHT-STER	All
456.499S	12MM/135 DEG TI CANN TROCH FIXATION NAIL 380MM/LEFT-STER	All
456.500S	12MM/135 DEG TI CANN TROCH FIXATION NAIL 400MM/RIGHT-STER	All
456.501S	12MM/135 DEG TI CANN TROCH FIXATION NAIL 400MM/LEFT-STER	All
456.502S	12MM/135 DEG TI CANN TROCH FIXATION NAIL 420MM/RIGHT-STER	All
456.503S	12MM/135 DEG TI CANN TROCH FIXATION NAIL 420MM/LEFT-STER	All
456.504S	12MM/135 DEG TI CANN TROCH FIXATION NAIL 440MM/RIGHT-STER	All
456.505S	12MM/135 DEG TI CANN TROCH FIXATION NAIL 440MM/LEFT-STER	All
456.506S	12MM/135 DEG TI CANN TROCH FIXATION NAIL 460MM/RIGHT-STER	All
456.507S	12MM/135 DEG TI CANN TROCH FIXATION NAIL 460MM/LEFT-STER	All
456.510S	12MM/125 DEG TI CANN TROCH FIXATION NAIL 235MM-STERILE	All
456.511S	12MM/130 DEG TI CANN TROCH FIXATION NAIL 235MM-STERILE	All
456.512S	12MM/135 DEG TI CANN TROCH FIXATION NAIL 235MM-STERILE	All
456.630S	14MM/130 DEG TI CANN TROCH FIXATION NAIL 300MM/RIGHT-STER	All
456.631S	14MM/130 DEG TI CANN TROCH FIXATION NAIL 300MM/LEFT-STER	All
456.632S	14MM/130 DEG TI CANN TROCH FIXATION NAIL 320MM/RIGHT-STER	All
456.633S	14MM/130 DEG TI CANN TROCH FIXATION NAIL 320MM/LEFT-STER	All
456.634S	14MM/130 DEG TI CANN TROCH FIXATION NAIL 340MM/RIGHT-STER	All
456.635S	14MM/130 DEG TI CANN TROCH FIXATION NAIL 340MM/LEFT-STER	All
456.636S	14MM/130 DEG TI CANN TROCH FIXATION NAIL 360MM/RIGHT-STER	All
456.637S	14MM/130 DEG TI CANN TROCH FIXATION NAIL 360MM/LEFT-STER	All
456.638S	14MM/130 DEG TI CANN TROCH FIXATION NAIL 380MM/RIGHT-STER	All
456.639S	14MM/130 DEG TI CANN TROCH FIXATION NAIL 380MM/LEFT-STER	All
456.640S	14MM/130 DEG TI CANN TROCH FIXATION NAIL 400MM/RIGHT-STER	All
456.641S	14MM/130 DEG TI CANN TROCH FIXATION NAIL 400MM/LEFT-STER	All

November 1st, 2012

Return Receipt Requested

SYNTHES GmbH
Luzernstrasse 21
CH-4528 Zuchwil
Tel +41 61 965 11 11

456.642S	14MM/130 DEG TI CANN TROCH FIXATION NAIL 420MM/RIGHT-STER	All
456.643S	14MM/130 DEG TI CANN TROCH FIXATION NAIL 420MM/LEFT-STER	All
456.644S	14MM/130 DEG TI CANN TROCH FIXATION NAIL 440MM/RIGHT-STER	All
456.645S	14MM/130 DEG TI CANN TROCH FIXATION NAIL 440MM/LEFT-STER	All
456.646S	14MM/130 DEG TI CANN TROCH FIXATION NAIL 460MM/RIGHT-STER	All
456.647S	14MM/130 DEG TI CANN TROCH FIXATION NAIL 460MM/LEFT-STER	All