

Datum: 2015-08-27
Utgåva: 1.0

Godkänd av: Jonas Holmgren
Sökväg:

Field Safety Notice

Manufacturer

Eurovema AB

Baldersvägen 38
SE-332 35 Gislaved
Sweden

Tel: +46-371-39 01 00
Fax: +46-371-189 82
E-post: info@eurovema.se

Contactperson:

Namn: Emma Nilsson

Tel: +46 734 42 59 32
E-post: emma.nilsson@eurovema.se

Product name or type: Miniflex, Flexmobil, Forma

Date: 2015-08-27

Actions to be taken: Replacement of component

Important!

Product reference:

See enclosed list

Problem or issue:

Listed power wheelchairs and work chairs are assembled with a seat cross containing an actuator bracket with manufacturing errors. A few components have broken.

Advice on actions to be taken by the distributor and the user:

Distributor: Examine your stock, if you find these products quarantine all stock. User of sold chairs must be contacted and date and place for replacement agreed. Contact manufacturer for replacement parts.

User: Determine your chair and if it is listed contact your national distributor.

Hereby we confirm that this issue has been reported to national competent authority.

Place/Date

Gislaved 2015-08-27



Emma Nilsson
Eurovema AB

Kvantitet / Qty	237
Qty Sweden	235
Qty Germany	2

Model	Serial number
--------------	----------------------

<i>Elrullstol Flexmobil</i>	1040098272
<i>Powe wheelchair Flexmobil</i>	1040101063
	1040121301
	1040121870
	1040132752
	1040135376

<i>Elrullstol Miniflex</i>	1040107201
<i>Power wheelchair Miniflex</i>	1040113362
	1040113473
	1040113825
	1040117329
	1040117347
	1040119056
	1040119142
	1040119415
	1040121863
	1040123012
	1040127091
	1040127203
	1040127205
	1040131361
	1040132735
	1040135490
	1040135658
	1040136996
	1040137334

<i>Arbetsstol Euroflex Forma EI</i>	104012533
<i>Work chair Euroflex Forma EI</i>	1040109002
	1040112015
	1040114359
	1040114360
	1040114641
	1040115993
	1040116960
	1040117002
	1040117321
	1040117323
	1040118121
	1040118488
	1040118491
	1040118493
	1040118683
	1040118809

1040119400
1040119403
1040119406
1040119409
1040119597
1040119600
1040119623
1040119871
1040119876
1040120480
1040120482
1040120568
1040121026
1040121114
1040121254
1040121864
1040122641
1040122729
1040125534
1040126189
1040126191
1040126499
1040126501
1040129975
1040130152
1040130622
1040130744
1040130885
1040131027
1040131558
1040131566
1040131945
1040132444
1040132739
1040132830
1040132832
1040132938
1040133970
1040134094
1040134126
1040134246
1040134247
1040134621
1040134959
1040135521
1040137021
1040137023
1040137025
1040137052
1040137064

1040137600
1040137605
1040137613
1040137615
1040137662
1040137664
1040122505
1040122532

Arbetsstol Euroflex Forma gas 1040113904
Work chair Euroflex Forma gas 1040114361
1040114362
1040114363
1040114364
1040115337
1040116885
1040116972
1040116978
1040116984
1040116991
1040117326
1040117364
1040117970
1040117971
1040117972
1040118229
1040118451
1040118765
1040118895
1040119183
1040119188
1040120111
1040120962
1040120963
1040120964
1040120965
1040120966
1040120997
1040120998
1040120999
1040121955
1040121956
1040121957
1040121958
1040122417
1040122426
1040122427
1040122530
1040122531
1040122643

1040122644
1040122645
1040123045
1040123297
1040123298
1040123807
1040125294
1040125311
1040125312
1040125313
1040125314
1040125315
1040126467
1040126605
1040126948
1040127045
1040127097
1040127362
1040127363
1040127364
1040127365
1040127366
1040128245
1040128420
1040128421
1040128422
1040129457
1040129733
1040129738
1040129743
1040130547
1040130899
1040130923
1040130944
1040131020
1040131212
1040131353
1040131538
1040131651
1040131669
1040132217
1040132222
1040132264
1040132265
1040132365
1040132378
1040132831
1040133345
1040133346
1040133347

1040133348
1040133349
1040133948
1040133976
1040133977
1040134111
1040134597
1040134616
1040134626
1040134861
1040134868
1040134946
1040135313
1040135314
1040135315
1040135316
1040135317
1040135355
1040135629
1040135636
1040136695
1040136696
1040136697
1040136698
1040136699
1040136711
1040136717
1040137534
1040137535
1040137536
1040137537
1040137538
1040137567
1040121959
1040123314
1040125519
104125525
1040129344
1040122666
1040125363
1040127925
1040126606
1040125525

Reservdel/Spare part

1040067141
1040128155